


NORMA Automated Cell Counters & Cell Viability Analyzers






 We are here

Discover our other range

CYTONOTE


Live Cell Imaging Inside the Incubator



IPRASENSE reinvents the automated Cell Counter and Viability Analyzers.

Our unique label-free imaging technology provides extremely fast cell count and viability from a few μ l sample volume of your cell suspension.

The unmatched repeatability directly results from the extremely large field of view of the single analyzed image, together with the sample preparation free method (no dilution, no label like trypan blue).



Several thousands of cells counted within a single image gives unmatched rapidity and repeatability



**SAMPLE
PREPARATION FREE**



**SHORT TIME
RESULTS**



**HIGH
REPETABILITY**



**LOW SAMPLE
VOLUME**

Mammalian Cell Culture - Cell Line - Media / Process Development - Drug Discovery



NORMA XS

Fast, simple & robust solutions for your routine lab cell counts.



NORMA 4S

Ready for fast, simple & 100% automatic solutions on your parallel bioreactors.



NORMA HT

Fast, simple & robust solutions for your routine lab cell counts and your high throughput platforms.

APPLICATIONS

- ✓ VIABILITY
- ✓ RATIO ASPECTS
- ✓ GROWTH CURVES
- ✓ MATCH WITH REFERENCE TRYPAN BLUE METHOD
- ✓ AUTOMATIC CELL COUNT
- ✓ CELL SIZE DISTRIBUTION

For research use only (RUO). Not for use in diagnostic procedures.

NORMA XS

Cell Counters & Viability Analyzers

The **NORMA XS** is the most simple automatic benchtop. It uses the new and revolutionary lensless imaging technology to measure total and viable cell count and determine viability. The **NORMA XS** provides an unmatched accuracy and precision thanks to its wide field of view that allows counting several thousands of cells within a single sample image.

The measurement results are returned by the **HORUS** software. Each parameter is visible on a simple graphical interface with the possibility to follow several running cultures on user friendly charts.




Image zoom of CHO cells with overlapped mask of cell detection and viability determination.

Green circles (Viable cells)
Red circles (Dead cells)

ABSOLUTE EASE OF USE




Take your unprepared sample.
(Label free, Undiluted)



Load the sample into the slide chamber.



Insert the slide into the NORMA XS.



Click "▶" in HORUS Software for launch measure.

- Cells** • Mammalian Cells
- Concentration range** • 10^4 - $4 \cdot 10^7$ cell/ml
- Cell size range** • 7 - 50 μ m
- Sample volume** • 3 - 13 μ l
- Numbers of sample** • 4 samples/slide - 100 samples/box
- Viability determination** • Light diffraction
- Counting time** • 10 to 15 seconds
- Image** • .PNG / .BMP / .TIFF / .RAW / .AVI
- Dimensions** • 15 x 7,5 x 10,5 cm
- Enclosure** • Stainless steel, POM,
- Weight** • 1 kg
- Power supply** • USB
- Pharmaceutical industries** • 21 CFR part 11 & IQ/OQ
- Integration** • No

CELL LINES EXPERIENCE WITH NORMA

MAMMALIAN CELLS

CHO	JURKAT
HEK 293	YT
NIH 3T3	PC-12
HELA	VERO

INSECT CELLS

SF9	HIGH FIVE
-----	-----------


NORMA 4S

Cell Counters & Viability Analyzers

THE NORMA 4S CELL COUNTER FITS PERFECTLY WITH APPLICATIONS REQUIRING AUTOMATION AND HIGH THROUGHPUT

The **NORMA 4S** is a fully automatic cell counter for high throughput cell culture monitoring. It rapidly measures cell concentration and viability without repeatability compromise. Each samples is analyzed undiluted.

The **NORMA 4S** is ready for receiving the sample from the robotic arms of a parallel micro bioreactor, an automatic sampler or even a manual pipette.



The **NORMA 4S** works with precise calibrated measurement chambers **1** constructed on single-use slides **2**. The refillable slide cartridge **3** is ready to run 144 samples **4** without user interaction.

Ambr15® integration



INTEGRATION IN HIGH THROUGHPUT
AUTOMATED MICRO BIOREACTOR

Cells	• Mammalian Cells
Concentration range	• 10 ⁴ -4.10 ⁷ cell/ml
Cell size range	• 7 - 50 µm
Sample volume	• 3 - 13 µl
Numbers of sample	• 4 samples/slide - 144 samples/box
Viability determination	• Light diffraction
Counting time	• 15 seconds
Image	• .PNG / .BMP / .TIFF / .RAW / .AVI
Dimensions	• 30 x 11,5 x 25 cm
Enclosure	• Stainless steel, POM,
Weight	• 4 kg
Power supply	• USB + 24 V DC (110 - 240 V AC)
Pharmaceutical industries	• 21 CFR part 11 & IQ/OQ
Integration	• Ambr15®

CELL LINES EXPERIENCE WITH NORMA

MAMMALIAN CELLS

CHO	JURKAT
HEK 293	YT
NIH 3T3	PC-12
HELA	VERO

INSECT CELLS


SF9	HIGH FIVE
-----	-----------

NORMA HT

Cell Counters & Viability Analyzers


THE NORMA HT CELL COUNTER IS DEDICATED TO CELL COUNT AND VIABILITY FOR MULTIWELL PLATES CULTURES

The **NORMA HT** cell counter is the most simple automatic benchtop cell counter for high throughput parallel culture monitoring. It rapidly measure viable cell count and viability on up to 24 samples. Each samples is analyzed undiluted and several thousands of cells are counted within a single image.





Each sample can be loaded manually to the 24 chambers counting slide with standard mono or multi channel pipettes or the complete sampling and counting procedure can be fully automated / integrated with liquid handlers and robotic systems.

Image zoom of CHO cells with overlapped mask of cell detection and viability determination.



Green circles (Viable cells)
Red circles (Dead cells)

HOW TO USE IT ?



The measurement results are returned by the **HORUS** software. Each parameter is visible on a simple graphical interface with the possibility to follow up to 24 running cultures on user friendly charts.

- Cells** • Mammalian Cells
- Concentration range** • 10^4 - $4 \cdot 10^7$ cell/ml
- Cell size range** • 7 - 50 μ m
- Sample volume** • 2 - 13 μ l
- Numbers of sample** • 4 samples/slide - 100 samples/box
- Viability determination** • Light diffraction
- Counting time** • 15 seconds
- Image** • .PNG / .BMP / .TIFF / .RAW / .AVI
- Dimensions** • 29,5 x 26,5 x 29,5 cm
- Enclosure** • Stainless steel,
- Weight** • 12 kg
- Power supply** • 100 - 240 V AC + USB
- Pharmaceutical industries** • 21 CFR part 11 & IQ/OQ
- Integration** • TECAN & Biomek

CELL LINES EXPERIENCE WITH NORMA

MAMMALIAN CELLS

CHO	JURKAT
HEK 293	YT
NIH 3T3	PC-12
HELA	VERO

INSECT CELLS

SF9	HIGH FIVE
-----	-----------

- Cells** • Mammalian Cells
- Concentration range** • 10^4 - $4 \cdot 10^7$ cell/ml
- Cell size range** • 7 - 50 μ m
- Sample volume** • 3 - 13 μ l
- Numbers of sample** • 4 samples/slide - 100 samples/box
- Viability determination** • Light diffraction
- Counting time** • 10 to 15 seconds
- Image** • .PNG / .BMP / .TIFF / .RAW / .AVI
- Dimensions** • 15 x 7,5 x 10,5 cm
- Enclosure** • Stainless steel, POM,
- Weight** • 1 kg
- Power supply** • USB
- Pharmaceutical industries** • 21 CFR part 11 & IQ/OQ
- Integration** • No

CELL LINES EXPERIENCE WITH NORMA

CHO	JURKAT
HEK 293	YT
NIH 3T3	PC-12
HELA	VERO
SF9	HIGH FIVE



- Cells** • Mammalian Cells
- Concentration range** • 10^4 - $4 \cdot 10^7$ cell/ml
- Cell size range** • 7 - 50 μ m
- Sample volume** • 3 - 13 μ l
- Numbers of sample** • 4 samples/slide - 144 samples/box
- Viability determination** • Light diffraction
- Counting time** • 15 seconds
- Image** • .PNG / .BMP / .TIFF / .RAW / .AVI
- Dimensions** • 30 x 11,5 x 25 cm
- Enclosure** • Stainless steel, POM,
- Weight** • 4 kg
- Power supply** • USB + 24 V DC (110 - 240 V AC)
- Pharmaceutical industries** • 21 CFR part 11 & IQ/OQ
- Integration** • Ambr15®

CELL LINES EXPERIENCE WITH NORMA

CHO	JURKAT
HEK 293	YT
NIH 3T3	PC-12
HELA	VERO
SF9	HIGH FIVE



- Cells** • Mammalian Cells
- Concentration range** • 10^4 - $4 \cdot 10^7$ cell/ml
- Cell size range** • 7 - 50 μ m
- Sample volume** • 2 - 13 μ l
- Numbers of sample** • 24 samples/slide - 480 samples/box
- Viability determination** • Light diffraction
- Counting time** • 15 seconds
- Image** • .PNG / .BMP / .TIFF / .RAW / .AVI
- Dimensions** • 29,5 x 26,5 x 29,5 cm
- Enclosure** • Stainless steel,
- Weight** • 12 kg
- Power supply** • 100 - 240 V AC + USB
- Pharmaceutical industries** • 21 CFR part 11 & IQ/OQ
- Integration** • TECAN & Biomek

CELL LINES EXPERIENCE WITH NORMA

CHO	JURKAT
HEK 293	YT
NIH 3T3	PC-12
HELA	VERO
SF9	HIGH FIVE



CONTACT

📍 5 Avenue de l'Europe,
Hélioparc,
34830 Clapiers,
France

☎ + 33 4 99 65 48 42

✉ info@iprasense.com

🌐 www.iprasense.com

